BONE SURGERY Ready To Use

en	Instructions for use	3
es	Instrucciones de uso	10
fr	Instructions d'utilisation	18
de	Gebrauchsanweisung	26

THERAVET

DESCRIPTION

BIOCERA-VET® BONE SURGERY Ready To Use (BIOCERA-VET® BS RTU) is a synthetic, biocompatible, osteoconductive and bioresorbable bone substitute material for filling non-load-bearing bone defects. The material is applied directly to the bone defect from the syringe or with a cannula without any further preparation.

BIOCERA-VET® BS RTU is a sterile single-use product, available in 2 different dosage forms (1 cc and 3 cc).

Thanks to its high injectability and self-setting (hardening), BIOCERA-VET® BS RTU is easy and intuitive to use.

COMPOSITION

BIOCERA-VET® BS RTU is a mineral bone cement composed of synthetic calcium and phosphate salts finely dispersed in a biocompatible carrier liquid made from short-chain triglycerides (caprylic/capric triglycerides) as well as two emulsifiers (polyoxyl 35 castor oil and cetyl phosphate). Caprylic/capric triglycerides and polyoxyl 35 castor oil are made from plant-based raw materials.

The setting reaction of BIOCERA-VET® BS RTU starts after application upon contact with an aqueous solution, such as body fluid. BIOCERA-VET® BS RTU sets in situ to a microcrystalline, calcium-deficient hydroxyapatite (CDHA) and alpha-tricalcium phosphate (α -TCP), which are the main phases. The minor phase consists of calcium hydrogen phosphate (monetite) and calcium carbonate (calcite). The chemical composition and crystalline structure of the final reaction product (CDHA) corresponds largely to the mineral components of natural bone.

BIOCERA-VET® BS RTU does not contain any substances of animal origin, added preservatives or pharmacologically active agents. A volume of 1 cc BIOCERA-VET® BS RTU contains approx. 24 mg potassium in the form of K₂HPO₄ (see dosage and application).

Table 1: Composition of BIOCERA-VET® BS RTU

Components	Percentage (%)
Alpha-tricalcium phosphate	48.35 - 49.88
Calcium hydrogen phosphate (monetite)	20.95 - 21.61
Calcium carbonate (calcite)	8.06 - 8.32
Tricalcium orthophosphate (CDHA)	3.23 - 3.33
Dipotassium hydrogen phosphate	2.41 - 2.49
Miglyol 812 (caprylic/capric triglycerides)	11.56 - 13.68
Kolliphor ELP (polyoxyl 35 castor oil)	2.11 - 2.50
Amphisol A (cetyl phosphate)	0.70 - 0.82

INDICATIONS

BIOCERA-VET® BS RTU is a self-setting, synthetic bone substitute material for filling non-infected bone defects in companion animals.

TARGET SPECIES

BIOCERA-VET® BS RTU is intended for bone void filling as add-on of bone surgeries in dogs and cats.

AREA OF APPLICATION

BIOCERA-VET® BS RTU is intended for filling or reconstructing non-load-bearing bone defects or for filling bone defects that have been sufficiently stabilized by appropriate means.

Particular areas of application are:

- Bone fractures with bone defects, e.g. fractures of the tibia, ulna, or femur
- bone fusion at a joint level – Arthrodesis e.g. carpal or tarsal joint
- bone defects following bone resection as corrective osteotomy in angular or rotational limb deformity
- bone defects following tibial tuberosity advancement (TTA)
- bone defects following resection of benign bone tumors or bone cysts
- bone defects after removal or replacement of osteosynthesis implants
- support for the fixation of osteosynthesis implants (e.g. of bone screws)

BIOCERA-VET® BS RTU must only be applied to a well vascularized and non-infected bone site.

The correct reduction and stabilization of fractures is to be assured by means of appropriate fixation.

INSTRUCTION FOR USE

BIOCERA-VET® BS RTU is an implantable product and designed for use in open or minimally invasive applications.

Fig. 1: BIOCERA-VET® BONE SURGERY
Ready To Use product presented in pre-filled syringe of 1 cc and 3 cc (ready to use)

The BIOCERA-VET® BS RTU package includes a cannula. When using the cannula, the force required for application is increased. It must be noted that, some of the BIOCERA-VET® BS RTU will remain in the cannula and is thus not available for filling the defect. BIOCERA-VET® BS RTU can also be applied without using a cannula.

The veterinarian is responsible for the patient's treatment plan, including the duration and timing of the clinical and radiological follow-up. The patient must follow to the veterinarian's treatment plan. During the pre-operation discussion, the patient's owner should be informed about the treatment conditions with BIOCERA-VET® BS RTU according to the instructions for use.

Preparation of the 1 cc and 3 cc version

The BIOCERA-VET® BS RTU package contains a pre-filled syringe and a cannula.

Open the package that contains the syringe with BIOCERA-VET®BS RTU and remove the blue cap. BIOCERA-VET® BS RTU can then be applied to the prepared bone defect directly or using the cannula. To do so, apply slow, and even pressure to the syringe plunger (see Figure 2 and 3).

The cannula can be shortened according to the requirements (e.g. with Cooper scissors, Figure 3).


Fig. 2: Application of cannula for 1cc Calcium Phosphate Cement Ready to Use


Fig. 3: Application of cannula for 3cc Calcium Phosphate Cement Ready to Use

Application

The amount of BIOCERA-VET® BS RTU required to fill a defect depends on the defect size. In order to avoid delays during the surgical procedure, it should be ensured prior to treatment that a sufficient number of BIOCERA-VET® BS RTU packages are available to fill the bone defect completely.

For minimally invasive applications, suitable imaging methods should be used to monitor the filling of the bone defect.

Setting behavior / hardening

BIOCERA-VET® BS RTU is set after combination or contact with an aqueous solution, such as body fluids. Prior to contact with an aqueous solution, BIOCERA-VET® BS RTU can be processed indefinitely.

Like all calcium phosphate cements, BIOCERA-VET® BS RTU must be injected only after final reduction and stabilization of the bone defect to prevent manipulation during the further hardening phase and impairing the setting behavior. Once it has been applied into the bone defect, BIOCERA-VET® BS RTU must not be manipulated (no dabbing, no corrective measures, no testing of the hardening, etc.). Since BIOCERA-VET® BS RTU is hardened with the surrounding fluid, the strength development depends on the shape and size of the filled defect.

The surface dimensional stability of BIOCERA-VET® BS RTU is achieved after approx. 10 to 15 minutes in an aqueous environment, after approx. 3 days, the compression strength reaches up to 35 MPa.

BIOCERA-VET® BS RTU takes several days to set in the implantation site and reaches a compressive strength of up to 45 MPa in a fully hardened state.

Consistency and compressive strength

After initiating the setting process the compressive strength of BIOCERA-VET® BS RTU increases steadily. Ten minutes after injection, a hard outer shell is formed covering the pasty core. 16 to 24 hours after injection, the compressive strength is comparable to cancellous bone. The ultimate compressive strength is three to four times higher than healthy cancellous bone and reached after about four days depending on the volume/geometry of the filled bone defect and body fluid supply.

DOSAGE

BIOCERA-VET® BS RTU is available in the following volumes:

Content per syringe		Potassium content per syringe
1cc	2g	Approx. 24 mg
3cc	6g	Approx. 70 mg

STERILIZATION

BIOCERA-VET® BS RTU is a sterile medical device. It is sterilized using gamma irradiation. Due to the risk of infection transmission and/or the potential impairment of product performance, BIOCERA-VET® BS RTU must not be cleaned or resterilized.

BIOCERA-VET® BS RTU is intended for single use only.

STORAGE AND SHELF LIFE

It is recommended to store BIOCERA-VET® BS RTU in a dry place at room temperature (between 5 to 25°C). The required application force and setting speed may change if stored outside of the stated temperature range.

Short-term deviations from the storage conditions during transportation (maximum of 5 days) have been tested for the range from – 23°C to +55°C and do not result in an unusable product.

As this is a sterile medical device, the packaging must be protected from damage to avoid contamination of the product itself.

Syringes from damaged packages must not be used. The expiry date is on the labels. The product must not be used after this date.

ADDITIONAL NOTIFICATIONS

Caution: Federal (USA) law restricts this device to be used by or on the order of a licensed veterinarian.

For animal use only.

DISPOSAL

No special disposal measures are required for unopened products. Explanted or contaminated material must be disposed of in accordance with the veterinary practices' standard disposal procedure.

CONTRAINDICATIONS

BIOCERA-VET® BS RTU must not be used in case of:

- acute or chronic infections at the implantation site, e.g. osteomyelitis
- bone defects adjacent to open epiphyseal plates
- cranioplasty
- pregnant or nursing animals
- known intolerance to an ingredient of BIOCERA-VET® BS RTU (see composition)

BIOCERA-VET® BS RTU must be used only after carefully weighing the risks and benefits in the case of:

- disorders of the bone metabolism
- endocrinopathies
- immunosuppressive therapy
- simultaneous therapy with medicines that influence bone metabolism

USE DURING PREGNANCY OR LACTATION

No clinical data for use in pregnant or lactating animal is available. Therefore, BIOCERA-VET® BS RTU should not be used during pregnancy or lactation.

USE IN PUPPIES

There are no clinical known reasons that would restrict use in puppies and young animals.

INTERACTIONS

If the patient is being treated simultaneously with resorption-inhibiting substances (particularly bisphosphonates, long term treatments with NSAIDs – non-steroidal anti-inflammatory drugs), slower resorption of the implant material should be assumed.

No additional interactions with other medical devices or medicinal products are known, as long as they do not have a direct effect on bone metabolism. BIOCERA-VET® BS RTU is a non-metallic, non-conducting and non-magnetic bone substitute material and therefore is designated as MRI safe.

PRECAUTIONS AND WARNINGS

The use of BIOCERA-VET® BS RTU is restricted to specialists who are familiar with handling bone substitute materials, the appropriate surgical techniques, and the treatment of bone defects from their training.

When using BIOCERA-VET® BS RTU, as with all surgical interventions, surgery-related risks may occur, particularly rejection reactions, pseudarthrosis, swelling, seroma and hematoma formation, pain, and wound healing disorders.

Subsequent revision operations may become necessary.

Patients with weaker immune systems are at a higher risk of infections and implant failure. Owners of such patients must be informed by medical staff of the possible risks before surgery.

The treatment of post-operative infections may be hampered by the presence of an implanted foreign body, and it may prove necessary to remove the implanted material.

Conventional surgical tools that are part of the standard equipment in the operating theatre can be used to remove bone cement.

Any implanted foreign body can lead to fever or an allergic or inflammation reactions.

In humans, very rare cases of allergic reactions and anaphylactic shock have been reported in literature associated with the component macrogolglycerol ricinoleate Ph. Eur. polyoxyl 35 castor oil, which is contained in BIOCERA-VET® BS RTU.

In dogs no toxic and allergic reactions have been reported after local injection of polyoxyl 35 castor oil.

For patients with severely impaired kidney or liver function or adrenal insufficiency, the additional intake of small amounts of potassium can increase the danger of hyperkalemia or aggravate the existing condition. This is also the case for patients with reduced renal potassium excretion induced by the use of medicines (e.g. heparin, ACE inhibitors, potassium-sparing diuretics, spironolactone, non-steroidal anti-inflammatory drugs, cyclosporin A).

Since potassium is only successively released from BIOCERA-VET® BS RTU and the amount of potassium contained in BIOCERA-VET® BS RTU corresponds to only a fraction of the amount present in the daily nutritional intake, a low risk can be assumed, even with severely impaired kidney function.

BIOCERA-VET® BS RTU must only be implanted after sufficient debridement of the bone defect to ensure a vital bone site. The defect must be completely filled in order to establish direct osseous contact between BIOCERA-VET® BS RTU and the surrounding bone.

When applying BIOCERA-VET® BS RTU into a bone defect, it must be ensured that the material does not escape into adjacent soft tissue or blood vessels. To prevent an embolism, it must be ensured that no material is inserted into open veins or arteries, particularly when inserting BIOCERA-VET® BS RTU under pressure into defects that are enclosed on all sides.

When BIOCERA-VET® BS RTU is applied into heavily bleeding bone defects, the material may be pushed out again through bleeding pressure. It is recommended to monitor bleeding before applying BIOCERA-VET® BS RTU and to continue filling the defect only when bleeding is under control.

Due to its mechanical properties, BIOCERA-VET® BS RTU can support the stabilization of bone defects, but the actual stabilization must be ensured by other means.

BIOCERA-VET® BS RTU must not be mixed with aqueous solutions, including those of autologous or allogenic origin (e.g. blood) as this results in changes to the material properties.

BIOCERA-VET® BS RTU is resorbed by biological processes and replaced by the body's own bone. BIOCERA-VET® BS RTU resorbs slowly. Depending on the implantation conditions and the metabolic activity at the implantation site, BIOCERA-VET® BS RTU can also remain permanently in the body as an osseous integrated material.

Any unused content of opened or damaged packages must not be used for further operations and must be discarded (see disposal).

SIDE EFFECTS

No side effects attributable to the application of BIOCERA-VET® BS RTU are known.

MANUFACTURING

Manufactured by:

INNOTERE

biomaterial

INNOTERE GmbH
Meissner Str. 191
01445 Radebeul
GERMANY
+49 351 2599 9410
www.innotere.de

For TheraVet SA:



TheraVet

TheraVet
32/1, Avenue Jean Mermoz
6041 GOSSELIES
BELGIUM
www.thera.vet

Customer Service

Email: customerservice@thera.vet
Phone: +32 71 96 00 44

SYMBOLS

	Single use
	Do not use if the sterile protection envelope is damaged
	Prior to use, read instructions for use
	Expiry date
	Sterile by irradiation
	Storage conditions
	Do not resterilize
	Read the user manual
	Double sterile barrier system
	MR safe
	Manufacturer
	Distributor

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