## OSTEOSARCOMA Ready to use

|    |                            |    |
|----|----------------------------|----|
| en | Instructions for use       | 2  |
| es | Instrucciones de uso       | 9  |
| fr | Instructions d'utilisation | 17 |
| de | Gebrauchsanweisung         | 25 |

# THERAVET

**DESCRIPTION**

**BIOCERA-VET® OSTEOSARCOMA Ready To Use (BIOCERA-VET® OSA RTU)** product is a synthetic, biocompatible, osteoconductive and bioresorbable bone substitute material for filling bone defect in weight-bearing bones. The material is applied directly into the bone defect from the syringe or with a canula without any further preparation.

BIOCERA-VET® OSA RTU is a sterile single-use product, available in two packaging sizes (12 cc and 2x 12 cc).

Thanks to its high injectability and self-setting (hardening), BIOCERA-VET® OSA RTU is easy and intuitive to use.

**COMPOSITION**

BIOCERA-VET® OSA RTU is a mineral bone cement composed of synthetic calcium and phosphate salts finely dispersed in a biocompatible carrier liquid made from short-chain triglycerides (caprylic/capric triglycerides) as well as two emulsifiers (polyoxyl 35 castor oil and cetyl phosphate). Caprylic/capric triglycerides and polyoxyl 35 castor oil are made from plant-based raw materials.

The setting reaction of BIOCERA-VET® OSA RTU starts after application upon contact with an aqueous solution, such as body fluid. BIOCERA-VET® OSA RTU sets in situ to a microcrystalline, calcium-deficient hydroxyapatite (CDHA) as the main phases. The minor phase consists of calcium hydrogen phosphate (monetite) and calcium carbonate (calcite). The chemical composition and crystalline structure of the final reaction product (CDHA) correspond largely to the mineral components of natural bone.

BIOCERA-VET® OSA RTU does not contain any substances of animal origin, added preservatives or pharmacologically active agents. A volume of 12 cc BIOCERA-VET® OSA RTU contains approx. 280 mg potassium in the form of K<sub>2</sub>HPO<sub>4</sub> (see dosage and application).

**Table 1: Composition of BIOCERA-VET® OSA RTU**

| Components                                  | Percentage (%) |
|---------------------------------------------|----------------|
| Alpha-tricalcium phosphate                  | 48.35 - 49.88  |
| Calcium hydrogen phosphate (monetite)       | 20.95 - 21.61  |
| Calcium carbonate (calcite)                 | 8.06 - 8.32    |
| Tricalcium orthophosphate (CDHA)            | 3.23 - 3.33    |
| Dipotassium hydrogen phosphate              | 2.41 - 2.49    |
| Miglyol 812 (caprylic/capric triglycerides) | 11.56 - 13.68  |
| Kolliphor ELP (polyoxyl 35 castor oil)      | 2.11 - 2.50    |
| Amphisol A (cetyl phosphate)                | 0.70 - 0.82    |

**INDICATIONS**

BIOCERA-VET® OSA RTU is recommended as a bone substitute for the cementoplasty of primary osteosarcoma tumors at the weight-bearing bones of the appendicular skeleton.

BIOCERA-VET® OSA RTU can also be used for other indications where the bone is weakened by a disease such as other bone tumors, metastasis, a traumatic or surgically induced injury, or for filling bone defect where large volume of bone substitute is required.

Cementoplasty is a low-traumatic and minimally invasive, surgical procedure performed in a short operating time that leads to a fast post operating recovery (Compared to the limb-sparing surgical procedure and the amputation).

Thanks to this surgical procedure, BIOCERA-VET® OSA RTU strengthens the injected bone thus enabling to preserve the function of the limb by clinically reducing lameness and the risk of pathological fractures. BIOCERA-VET® OSA RTU also helps to relieve pain, therefore increasing the animal's well-being by mechanically reinforcing the bone affected by the tumor.

### **TARGET SPECIES**

BIOCERA-VET® OSA RTU is intended for the cementoplasty at weightbearing bones of the appendicular skeleton of dogs and cats.

### **AREA OF APPLICATION**

BIOCERA-VET® OSA RTU is intended for cementoplasty of bone weakened by a disease, a bone injury (traumatic or surgically induced) and for filling large bone defect.

BIOCERA-VET® OSA RTU must only be applied to a well vascularized and non-infected bone site.

### **INSTRUCTION FOR USE**

BIOCERA-VET® OSA RTU is an implantable product and designed for use in open or minimally invasive applications.

**Fig. 1:** BIOCERA-VET® OSA RTU product presented in pre-filled syringe of 12 cc


The BIOCERA-VET® OSA RTU package includes a 7G cannula. When using the cannula, the force required for application is increased. It must be noted that, some of the BIOCERA-VET® OSA RTU will remain in the cannula and is thus not available for filling the defect.

BIOCERA-VET® OSA RTU can also be applied without using a cannula.

The veterinarian is responsible for the patient's treatment plan, including the duration and timing of the clinical and radiological follow-up. The patient must follow to the veterinarian's treatment plan. During the pre-operation discussion, the patient's owner should be informed about the treatment conditions with BIOCERA-VET® OSA RTU according to the instructions for use.


### **Preparation of 12 cc version**

The BIOCERA-VET® OSA RTU package contains a pre-filled syringe, a spindle drive and a cannula. Open the packages that contain the syringe with BIOCERA-VET® OSA RTU..

**Fig. 2:** Preparation of the syringe with spindle drive

Then open the package with the spindle drive (spindle nut and threaded spindle). The spindle nut must be pushed onto the back end of the barrel until a clear click sound is heard. Check that both sides of the spindle nut are locked. Screw the threaded spindle in the spindle nut until it is in contact with the piston. Remove the blue cap from the syringe. BIOCERA-VET® OSA RTU can then be applied to the prepared bone defect directly or using the cannula by slowly rotating the spindle clockwise as presented below.

**Fig. 3:** Application of the cannula (optional) and the 12cc syringe

### **Application**

The amount of BIOCERA-VET® OSA RTU required to fill a defect depends on the defect size. In order to avoid delays during the surgical procedure, it should be ensured prior to treatment that a sufficient number of BIOCERA-VET® OSA RTU packages are available to fill the bone defect completely. For minimally invasive applications, suitable imaging methods should be used to monitor the filling of the bone defect.

### **Setting behavior / hardening**

BIOCERA-VET® OSA RTU is set after combination or contact with an aqueous solution, such as body fluids. Prior to contact with an aqueous solution, BIOCERA-VET® OSA RTU can be processed indefinitely.

Like all calcium phosphate cements, BIOCERA-VET® OSA RTU must be injected only after final reduction and stabilization of the bone defect (if applicable) to prevent manipulation during the further hardening phase and impairing the setting behavior. Once it has been injected into the bone cavity, BIOCERA-VET® OSA RTU must not be manipulated (no dabbing, no corrective measures, no testing of the hardening, etc.).

Since BIOCERA-VET® OSA RTU is hardened with the surrounding fluid, the strength development depends on the shape and size of the filled defect. The surface dimensional stability of BIOCERA-VET® OSA RTU is achieved after approx. 10 to 15 minutes in an aqueous environment, after approx. 3 days, the compression strength reaches up to 35 MPa.

BIOCERA-VET® OSA RTU takes several days to set in the implantation site and reaches a compressive strength of up to 45 MPa in a fully hardened state.

### **Consistency and compressive strength**

After initiating the setting process the compressive strength of BIOCERA-VET® OSA RTU increases steadily. Ten minutes after injection, a hard outer shell is formed covering the pasty core. 16 to 24 hours after injection, the compressive strength is comparable to cancellous bone. The ultimate compressive strength is three to four times higher than healthy cancellous bone and reached after about four days depending on the volume/geometry of the filled bone defect and body fluid supply.

### **DOSAGE**

BIOCERA-VET® OSA RTU is available in the following volumes:

| Content per syringe |      | Potassium content |
|---------------------|------|-------------------|
| 12 cc               | 24 g | Approx. 280 mg    |
| 2x 12 cc            | 48 g | Approx. 560 mg    |

### **STERILIZATION**

BIOCERA-VET® OSA RTU is a sterile medical device. It is sterilized using gamma irradiation. Due to the risk of infection transmission and/or the potential impairment of product performance BIOCERA-VET® OSA RTU must not be cleaned or resterilized. BIOCERA-VET® OSA RTU is intended for single use only.

### **STORAGE AND SHELF LIFE**

It is recommended to store BIOCERA-VET® OSA RTU in a dry place at room temperature (between 5 to 25°C). The required application force and setting speed may change if stored outside of the stated temperature range. Short-term deviations from the storage conditions during transportation (maximum of 5 days) have been tested for the range from -23°C to +55°C and do not result in an unusable product.

As this is a sterile medical device, the packaging must be protected from damage to avoid contamination of the product itself. Syringes from damaged packages must not be used.

The expiry date is on the labels. The product must not be used after this date.

## **ADDITIONAL NOTIFICATIONS**

Caution: Federal (USA) law restricts this device to be used by or on the order of a licensed veterinarian.

For animal use only.

## **DISPOSAL**

No special disposal measures are required for unopened products. Explanted or contaminated material must be disposed of in accordance with the veterinary practices' standard disposal procedure.

## **CONTRAINDICATIONS**

BIOCERA-VET® OSA RTU must not be used in case of:

- acute or chronic infections at the implantation site, e.g. osteomyelitis
- bone defects adjacent to open epiphyseal plates
- cranioplasty
- pregnant or nursing animals
- known intolerance to an ingredient of BIOCERA-VET® OSA RTU (see composition)

BIOCERA-VET® OSA RTU must be used only after carefully weighing the risks and benefits in the case of:

- disorders of the bone metabolism
- endocrinopathies
- immunosuppressive therapy
- simultaneous therapy with medicines that influence bone metabolism

## **USE DURING PREGNANCY OR LACTATION**

No clinical data for use in pregnant or lactating animal is available. Therefore, BIOCERA-VET® OSA RTU should not be used during pregnancy or lactation.

## **USE IN PUPPIES**

There are no clinical known reasons that would restrict use in puppies and young animal.

## **INTERACTIONS**

If the patient is being treated simultaneously with resorption-inhibiting substances (particularly bisphosphonates, long term treatments with NSAIDs – non-steroidal anti-inflammatory drugs), slower resorption of the implant material should be assumed. No additional interactions with other medical devices or medicinal products are known, as long as they do not have a direct effect on bone metabolism.

BIOCERA-VET® OSA RTU is a non-metallic, non-conducting and non-magnetic bone substitute material and therefore is designated as MRI safe.

## **PRECAUTIONS AND WARNINGS**

The use of BIOCERA-VET® OSA RTU is restricted to specialists who are familiar with handling bone substitute materials, the appropriate surgical techniques, and the treatment of bone defects from their training.

When using BIOCERA-VET® OSA RTU, as with all surgical interventions, surgery-related risks may occur, particularly rejection reactions, pseudarthrosis, swelling, seroma and

hematoma formation, pain, and wound healing disorders. Subsequent revision operations may become necessary.

Patients with weaker immune systems are at a higher risk of infections and implant failure. Owner of such patients must be informed by medical staff of the possible risks before surgery.

The treatment of post-operative infections may be hampered by the presence of an implanted foreign body and it may prove necessary to remove the implanted material.

Conventional surgical tools that are part of the standard equipment in the operating theatre can be used to remove bone substitute.

Any implanted foreign body can lead to fever or an allergic or inflammation reactions.

In humans, very rare cases of allergic reactions and anaphylactic shock have been reported in literature associated with the component macrogolglycerol ricinoleate Ph. Eur. polyoxyl 35 castor oil, which is contained in BIOCERA-VET® OSA RTU. In dogs, no toxic and allergic reactions have been reported after local injection of polyoxyl 35 castor oil.

For patients with severely impaired kidney or liver function or adrenal insufficiency, the additional intake of small amounts of potassium can increase the danger of hyperkalemia or aggravate the existing condition. This is also the case for patients with reduced renal potassium excretion induced by the use of medicines (e.g. heparin, ACE inhibitors, potassium-sparing diuretics, spironolactone, non-steroidal anti-inflammatory drugs, cyclosporin A). Since potassium is only successively released from BIOCERA-VET® OSA RTU and the amount of potassium contained in BIOCERA-VET® OSA RTU corresponds to only a fraction of the amount present in the daily nutritional intake, a low risk can be assumed, even with severely impaired kidney function.

BIOCERA-VET® OSA RTU must only be implanted after sufficient debridement of the bone defect in order to ensure a vital bone site. The defect must be completely filled in order to establish direct osseous contact between BIOCERA-VET® OSA RTU and the surrounding bone.

When injecting BIOCERA-VET® OSA RTU into a bone defect, it must be ensured that the material does not escape into adjacent soft tissue or blood vessels. To prevent an embolism, it must be ensured that no material is inserted into open veins or arteries, particularly when inserting BIOCERA-VET® OSA RTU under pressure into defects that are enclosed on all sides.

When BIOCERA-VET® OSA RTU is injected into heavily bleeding bone defect, the material may be pushed out again through bleeding pressure. It is recommended to monitor bleeding before applying BIOCERA-VET® OSA RTU and to continue filling the defect only when bleeding is under control.

Due to its mechanical properties, BIOCERA-VET® OSA RTU can support the stabilization of bone, but the actual stabilization must be ensured by other means.

BIOCERA-VET® OSA RTU must not be mixed with aqueous solutions, including those of autologous or allogenic origin (e.g. blood) as this results in changes to the material properties.

BIOCERA-VET® OSA RTU is resorbed by biological processes and replaced by the body's own bone. BIOCERA-VET® OSA RTU resorbs slowly. Depending on the implantation conditions and the metabolic activity at the implantation site,

BIOCERA-VET® OSA RTU can also remain permanently in the body as an osseous integrated material.  
 Any unused content of opened or damaged packages must not be used for further operations and must be discarded (see disposal).

### SIDE EFFECTS

No side effects attributable to the application of BIOCERA-VET® OSA RTU are known.

### MANUFACTURING

Manufactured by:

#### INNOTERE

*biomaterial*

INNOTERE GmbH  
 Meissner Str. 191  
 01445 Radebeul  
 GERMANY  
 +49 351 2599 9410  
 www.innotere.de

For TheraVet SA:




#### TheraVet

TheraVet  
 32/1, Avenue Jean Mermoz  
 6041 GOSSELIES  
 BELGIUM  
 www.thera.vet

### Customer Service

Email: customerservice@thera.vet  
 Phone: +32 71 96 00 44

### SYMBOLS

|                                                                                     |                                                          |
|-------------------------------------------------------------------------------------|----------------------------------------------------------|
|  | Single use                                               |
|  | Do not use if the sterile protection envelope is damaged |
|  | Prior to use, read instructions for use                  |
|  | Expiry date                                              |
|  | Sterile by irradiation                                   |
|  | Storage conditions                                       |
|  | Do not resterilize                                       |
|  | Read the user manual                                     |
|  | Double sterile barrier system                            |
|  | MR safe                                                  |
|  | Manufacturer                                             |
|  | Distributor                                              |

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