## GRANULES

en	Instructions for use	2
es	Instrucciones de uso	7
fr	Instructions d'utilisation	12
de	Gebrauchsanweisung	17

# THERAVET

### DESCRIPTION

**BIOCERA-VET® GRANULES** is a synthetic, porous, biocompatible and bioresorbable bone substitute material for filling non-load-bearing bone defects.

BIOCERA-VET® GRANULES is a sterile single-use product.

BIOCERA-VET® GRANULES is available in different volumes and is presented as multipacks.

### COMPOSITION

BIOCERA-VET® GRANULES is a porous, mineral bone substitute material of synthetic calcium and phosphates. The calcium phosphate phases are present with microcrystalline, calcium-deficient hydroxyapatite (CDHA) and alpha-tricalcium phosphate ( $\alpha$ -TCP), which are the main phases. The minor phase consists of calcium hydrogen phosphate (monetite) and calcium carbonate (calcite). The chemical composition and crystalline structure of the final reaction product (CDHA) corresponds largely to the mineral components of natural bone.

### INDICATIONS

BIOCERA-VET® GRANULES is a synthetic, porous bone substitute material for filling non-infected bone defects in companion animals.

### TARGET SPECIES

BIOCERA-VET® GRANULES is intended for bone void filling as add-on of bone surgeries in dogs and cats.

### AREA OF APPLICATION

BIOCERA-VET® GRANULES is intended for filling or reconstructing non-load-bearing bone defects or for filling bone defects that have been sufficiently stabilized by appropriate means.

Particular areas of application are:

- bone fractures with bone defect e.g. fractures of the tibia, ulna, or femur
- bone fusion at a joint level – Arthrodesis e.g. carpal or tarsal joint
- bone defects following bone resection as corrective osteotomy in angular or rotational limb deformity
- bone defect following tibial tuberosity advancement (TTA)
- bone defects following resection of benign bone tumors or bone cysts
- filling of bone defects in oral and maxillofacial surgery
- filling of spinal cages

BIOCERA-VET® GRANULES must only be applied into a well vascularized and non-infected bone bed, in order to ensure the formation of new bone. The correct reduction and stabilization of fractures is to be assured through sufficient fixation.

BIOCERA-VET® GRANULES is a single-use implant material.

### INSTRUCTION FOR USE

BIOCERA-VET® GRANULES is an implantable product and designed for use in open surgical applications.

The BIOCERA-VET® GRANULES package contains the sterile, ready-to-use product.

BIOCERA-VET® GRANULES can be combined intra-operatively with allogeneic or autologous materials, particularly blood, blood-based products, bone marrow aspirate or autologous cancellous bone. In these cases, particular attention should be paid to ensuring aseptic conditions.

BIOCERA-VET® GRANULES must only be implanted after sufficient debridement of the bone defect in order to provide a vital bone site. The defect should be filled completely in order to achieve a direct osseous contact between BIOCERA-VET® GRANULES and the surrounding bone. If it is impossible to fill the defect completely with BIOCERA-VET® GRANULES alone, the remaining defect sites could be filled with autologous bone or allogeneic materials.

BIOCERA-VET® GRANULES is resorbed and replaced by the body's own tissues. Residual material may still be visible on X-rays for up to 12 months. Incomplete resorption does not constitute an impairment of its function as a defect filler. The speed of new bone formation is not affected.

## **DOSAGE AND APPLICATION**

BIOCERA-VET® GRANULES is available in the following granules size and volumes:

<b>Granules size</b>	<b>Content</b>
500-710 µm	5 x 1cc
	3 x 3cc
	3 x 5cc

No upper dose limit for the implantation of BIOCERA-VET® GRANULES can be fixed. The amount of BIOCERA-VET® GRANULES necessary to fill a defect is dependent on the defect size.

## **STERILIZATION**

BIOCERA-VET® GRANULES is a sterile medical device. Sterilization is achieved by means of gamma irradiation. Due to the risk of infection transmission and/or the potential impairment of product performance, BIOCERA-VET® GRANULES must not be cleaned or resterilized.

BIOCERA-VET® GRANULES is intended for single use only.

## **STORAGE AND SHELF LIFE**

No particular storage conditions must be observed. It is recommended to store BIOCERA-VET GRANULES in a dry place at room temperature. As it is a sterile medical device, the package has to be protected from damage in order to avoid contamination of the product itself.

BIOCERA-VET® GRANULES from damaged packages must not be used. The expiry date is on the labels. The product must not be used after this date.

## **ADDITIONAL NOTIFICATIONS**

Caution: Federal (USA) law restricts this device to be used by or on the order of a licensed veterinarian.

For animal use only.

## **DISPOSAL**

No special disposal measures are required for unopened products. Explanted or contaminated material must be disposed of in accordance with the veterinary practices' standard disposal procedure.

## **CONTRAINDICATIONS**

BIOCERA-VET® GRANULES must not be used in case of:

- acute or chronic infections at the implantation site, e.g. osteomyelitis
- bone defects in the area of open epiphyseal plates
- patients with known calcium metabolism disorders, e.g. hypercalcemia

BIOCERA-VET® GRANULES must be used only after carefully weighing the risks and benefits in the case of:

- disorders of the bone metabolism
- endocrinopathies
- immunosuppressive therapy
- simultaneous therapy with medicines that influence bone metabolism
- simultaneous treatment with resorption-inhibiting active ingredients (particularly bisphosphonates, long term treatments with NSAIDs – non-steroidal anti-inflammatory drugs)

## **USE DURING PREGNANCY OR LACTATION**

Studies in pregnant or nursing animals have not been conducted.

## **USE IN PUPPIES**

Studies in puppies have not been conducted. There are no known reasons which would restrict use in young animals.

## **INTERACTIONS**

If the patient is being treated at the same time with resorption-inhibiting active substances (particularly bisphosphonates, NSAIDs – non-steroidal anti-inflammatory drugs), it should be assumed that the resorption of the implant material will be slowed, an increased risk of subsequent bone breakage cannot be excluded.

No additional interactions with other medical devices or medicines are known, as long as they do not have a direct effect on bone metabolism (see Contraindications).

BIOCERA-VET® GRANULES is a non-metallic, non-conducting and non-magnetic scaffold and therefore is designated as MRI safe.

## **PRECAUTIONS AND WARNINGS**

The use of BIOCERA-VET® GRANULES is restricted to specialists who through their training are familiar with handling bone substitute materials, with the appropriate surgical techniques as well as with the treatment of bone defects.

Upon application of BIOCERA-VET® GRANULES, as with all surgical interventions, surgery-related risks may occur, particularly rejection reactions, pseudarthrosis, swelling, seroma and hematoma formation, pain, and wound healing disorders. Subsequent revision operations may become necessary.

The treatment of post-operative infections may be hampered by the presence of an implanted foreign body and it may prove necessary to remove the implanted material.

Any implanted foreign body can lead to fever or an allergic reactions or inflammation.

Very rare cases of allergic reactions and anaphylactic shock to the component Macrogolglycerol ricinoleate Ph. Eur., which is used in the manufacturing of BIOCERA-VET® GRANULES, are described in the literature. Despite the steps to remove this during the manufacturing process, traces of this substance may still be present in the product.

In order to avoid delays during the surgical procedure, it should be ensured prior to treatment that a sufficient number of BIOCERA-VET® GRANULES units is available, e.g. in the form of multipacks, to fill the bone defect completely. Any unused content of opened or damaged packages must not be used for further operations and must be discarded.

For patients with severely impaired kidney function, adrenal insufficiency or liver cirrhosis, the additional intake of small amounts of potassium can increase the danger of hyperkalemia or aggravate the existing condition. This is also the case for patients with reduced renal potassium excretion induced by the use of medicinal products (e.g. heparin, ACE inhibitors, potassium-sparing diuretics, spironolactone, non-steroidal anti-inflammatory drugs, cyclosporin A). In particular when large amounts of implanted BIOCERA-VET® GRANULES are used, serum potassium levels must be monitored.

Patients with weaker immune systems are at higher risks of infections and implant loss. Owner of such patients should be informed by medical staff of the possible risks before surgery.

## **SIDE EFFECTS**

No serious side effects attributable to BIOCERA-VET® GRANULES are known.

## **MANUFACTURING**

**Manufactured by:**

**INNOTERE**

*biomaterial*

INNOTERE GmbH  
Meissner Str. 191  
01445 Radebeul  
GERMANY  
+49 351 2599 9410  
www.innotere.de

**For TheraVet SA:**

**THERAVET**


TheraVet  
32/1, Avenue Jean Mermoz  
6041 GOSSELIES  
BELGIUM  
www.thera.vet

## Customer Service

Email: [customerservice@thera.vet](mailto:customerservice@thera.vet)

Phone: +32 71 96 00 44

## SYMBOLS

	Single use
	Do not use if the sterile protection envelope is damaged
	Prior to use, read instructions for use
	Expiry date
	Sterile by irradiation
	Do not resterilize
	Read the user manual
	Double sterile barrier system
	MR safe
	Manufacturer
	Distributor

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