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This report provides essential information concerning safety and clinical performance of CP-Scaffold. Please note that the following information is intended for patients and lay persons. It does not replace the instructions for use or the implantation card. Further, this report does not offer general advice for treating medical conditions. If you have any questions about your medical condition or the use of CP-Scaffold or if you believe you are experiencing any side effects associated with CP-Scaffold, we strongly recommend contacting your healthcare professional. Medical terms have been avoided for easier understanding or are mentioned in brackets.

## 1. General information on the product

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### Year of market launch (CE-certificate)

2017

## 2. Application of CP-Scaffold

CP-Scaffold is a scaffold used to fill bone defects.

### 2.1 Medical treatment

CP-Scaffold is intended for filling bone defects in non-load-bearing bone. If CP-Scaffold is to be used in load-bearing bone, the bone must be sufficiently stabilized by suitable fixation.

Fields of application for CP-Scaffold are:

- Fractures of long bones (e.g. fractures of the radius, tibia or humerus)
- Corrections of bone malposition (osteotomy)
- Bone defects after removal or replacement of implants

### 2.2 Patient population

CP-Scaffold is suitable for adults.

### 2.3 Contraindications

CP-Scaffold must not be used on:

- Infections of the bone at the implant site (osteomyelitis)
- Bone defects due to malignant tumours
- Bone defects in the area between the shaft and the end of a bone (epiphyseal joint)
- Disorders of calcium metabolism (e.g. hypercalcemia)
- Pregnant and breastfeeding women

The physician critically assesses the use of CP-Scaffold, if the following diseases are present:

- Disorders of growth processes in the bone (bone metabolism)
- Diseases of the endocrine system (endocrinopathy)
- Taking medications for the treatment of autoimmune diseases (immunosuppressants)
- Taking medications that affect bone metabolism
- If the patient is a child

## 3. Product description

CP-Scaffold is a porous shaped body composed mainly of calcium phosphates (see Table 1). Since hydroxyapatite is the solid component of bone, CP-Scaffold is easily recognized by and incorporated into bone. Over time, CP-Scaffold is broken down by the body and new bone is formed at the same time. The time it takes for CP-Scaffold to completely degrade depends on various factors, e.g. the size of the inserted scaffold as well as the age of the patient.

Table 1: Components of CP-Scaffold

| Components                              |       |
|-----------------------------------------|-------|
| calcium-deficient hydroxyapatite (CDHA) | ≥ 75% |
| alpha-tricalcium phosphate (α-TCP)      |       |
| calcium hydrogen phosphate (monetite)   | ≤ 25% |
| calcium carbonate (calcite)             |       |

CP-Scaffold is available in different shapes (block, cylinder and wedge) and sizes (see Figure 1). The shape and amount of CP-Scaffold required to fill a defect depends on the type and size of the bone defect.


Figure 1: Shapes of CP-Scaffold (block, cylinder, wedge)

CP-Scaffold is clearly visible on the X-ray after implantation. But as soon as the degradation process of CP-Scaffold begins, the visibility decreases.

CP-Scaffold does not contain any pharmaceuticals and no substances of human or animal origin.

CP-Scaffold is intended for single use on a single person.

#### 4. Risks and warnings

Please contact your doctor if you experience any side effects related to CP-Scaffold or if you are concerned about potential risks. This report is not a substitute for talking to your doctor.

##### 4.1 Residual risks and side effects

Common side effects of surgery include bleeding, bruising (hematomas), accumulation of fluids (seromas), swelling, fever, wound healing problems, infection, delayed or abnormal fracture healing (pseudarthrosis), and pain.

Allergic reactions or rejections may also occur with the use of bone graft substitutes. These reactions have never been observed for CP-Scaffold. With the wedge-shaped products, fracture has been seen on X-ray in rare cases. However, fracture of the bone substitute material does not pose a clinical risk, as the material does not perform a load-bearing function and the bone is stabilized by other measures.

##### 4.2 Warnings and precautions

The physician has the duty to inform you about possible risks before the operation. Below you can read the warnings and precautions in the instructions for use of CP-Scaffold. Please consult your physician if you have any questions.

- Patients with a weakened immune system (e.g. rheumatics or diabetics) as well as smokers and alcohol abusers are at increased risk of infections and implant failure.
- CP-Scaffold may contain very small amounts of a castor oil (polyoxyl-35 castor oil) for which very rare cases of allergic reactions and hypersensitivity reaction of the immune system (anaphylactic shock) have been described in the literature.
- Treatment of infections after surgery may require additional surgery to remove CP-Scaffold from your body.
- The surgeon fills your bone defect completely with CP-Scaffold. If complete filling of the defect is not possible, the remaining defect can be filled with the body's own bone or suitable foreign bone.
- CP-Scaffold is a bone substitute and can only support the stability of the bone. The bone defect itself is stabilized by other implants (e.g. plates or screws).
- In rare cases, CP-Scaffold may fracture in the bone defect. The fracture is visible on the X-ray image. The healing process of the bone is not negatively influenced by this.

- CP-Scaffold can be mixed with the body's own material during surgery, especially with blood, blood products, cells from the bone marrow (bone marrow aspirate) or through the spongy inner tissue of the bone (cancellous bone).
- CP-Scaffold is degraded (dissolved) by biological processes and replaced by your own bone. The duration of the degradation process depends on many factors, e.g. the size of CP-Scaffold and the condition of your bone. In rare cases, CP-Scaffold may also remain permanently in your bone.

## 5. Clinical evaluation

CP-Scaffold belongs to a group of bone graft substitutes that have been successfully used in clinical applications for decades. CP-Scaffold itself has been in clinical use since 2017. The preferred use of calcium phosphates to produce bone substitute materials is derived from their close similarity to the mineral component of bone.

INNOTERE has conducted a clinical study on the use of CP-Scaffold (wedge shape) in a correction of bone malposition (realignment osteotomy). In this study, the bone healing and resorption of CP-Scaffold was investigated in comparison to the product, OSferion. OSferion also belongs to the group of calcium phosphate-containing bone graft substitutes and has been on the market since 1999. The aim of the study was to prove that CP-Scaffold is at least as good as the established product OSferion. The study included 71 patients. The statistically validated results of the study demonstrate that the bone healing of CP-Scaffold is not inferior to the comparator product, OSferion. No failure of bone healing was observed in any patient. CP-Scaffold was well incorporated into the bone in every case. INNOTERE 3D Scaffold is also non-inferior to the comparator product, OSferion, in terms of resorption. Four adverse events (pain, fracture of the bone substitute, hinge fracture, wound healing disorder) and one infection were observed during the study. Except for the fracture of the bone substitute, these adverse events are not related to the use of CP-Scaffold but are due to general risks of surgery.

A comprehensive clinical evaluation report is available for CP-Scaffold, demonstrating the safety and clinical performance of CP-Scaffold. The evidence for this assessment report comes from literature databases and the clinical study described, as well as feedback from clinical practice. There is sufficient evidence of clinical benefit. No unacceptable, undesirable side effects, as well as residual risks, were identified. The benefit-risk-ratio can be clearly assessed for the application of the CP-Scaffold material.

## 6. Therapeutic alternatives

CP-Scaffold is one of the synthetic bone substitute materials for which the following therapeutic alternatives are available:

- the use of the patient's own bone (autologous bone replacement)
- the use of human foreign bone (allogenic bone substitute)
- the use of other synthetic bone substitute materials made of hydroxyapatite or other compositions (e.g. calcium sulphates)
- the combination of the above-mentioned bone substitute materials
- omission of a bone substitute material (if necessary for smaller defects)

The healing of a bone defect is influenced by many factors, such as the defect size, the anatomical location of the bone defect, the patient's health condition and underlying diseases. If you are considering alternative treatment methods, please consult your physician, as he or she can best consider your individual situation.

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